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---

### **Literatur**

1. Loveday H.P, Wilson J A, Pratta R J, Golsorkhia M, Tinglea A, Baka A, Brownea J, Prietob J, Wilcox M. epic3: „National Evidence-Based Guidelines for Preventing Healthcare- Associated Infections in NHS Hospitals in England“. Journal of Hospital Infection 86S1 (2014) S1–S70
2. Filon FL, Cerchi R. „Epidemiology of latex allergy in healthcare workers.“ Med Lav. 2008;99(2):108-112
3. Sussman G, Gold M. „Guidelines for the management of latex allergies and safe latex use in health care facilities“. <http://www.aaaai.org/allergist/allergies/Types/latex-allergy/Pages/latex-allergies-safe-use.aspx>. Zuletzt aufgerufen am 11. Dezember 2013.

Ansell,\* and <sup>™</sup> are trademarks owned by Ansell Limited or one of its affiliates. © 2017 All Rights Reserved.

#### **North America**

Ansell Healthcare Products LLC  
111 Wood Avenue South  
Suite 210  
Iselin, NJ 08830, USA

#### **Europe, Middle East & Africa**

Ansell Healthcare Europe NV  
Riverside Business Part  
Blvd International 55  
1070 Brussels, Belgium

#### **Asia Pacific**

Ansell Services Asia Sdn. Bhd.  
Prima 6, Prima Avenue  
Block 3512, Jalan Teknokrat 6  
63000 Cyberjaya, Malaysia

#### **Australia & New Zealand**

Ansell Limited  
Level 3, 678 Victoria Street  
Richmond, Vic, 3121  
Australia

4. Charous BL, Hamilton RG, Yunginger JW. „Occupational latex exposure: characteristics of contact and systemic reactions in 47 workers.“ J Allergy Clin Immunol. 1994;94(1):12-18.
5. FDA Centre for Devices and Radiological Health. „Guidance for Industry and FDA Reviewers/Staff Premarket Notification“ [510(k)] „Submissions for Testing for Skin Sensitization to Chemicals in Natural Rubber Products“, Januar 1999
6. Geier J., Lessmann H., Mahler V., Pohrt U., Uter W. und Schnuch A., „Occupational contact allergy caused by rubber gloves - nothing has changed“. Contact Dermatitis, 67, 149–156 (2012)
7. Lucas D. Ansell Research and Development Melaka, Malaysia. 28. April 2016
8. Taylor P. „Chemical Allergy in Healthcare“. Asian Hospital & Healthcare Management. <http://www.asianhnm.com/medical-sciences/chemical-allergy>.  
Letzter Zugriff am 22. April 2016.

Ansell,® and <sup>™</sup> are trademarks owned by Ansell Limited or one of its affiliates. © 2017 All Rights Reserved.

**North America**

Ansell Healthcare Products LCC  
111 Wood Avenue South  
Suite 210  
Iselin, NJ 08830, USA

**Europe, Middle East & Africa**

Ansell Healthcare Europe NV  
Riverside Business Part  
Blvd International 55  
1070 Brussels, Belgium

**Asia Pacific**

Ansell Services Asia Sdn. Bhd.  
Prima 6, Prima Avenue  
Block 3512, Jalan Teknokrat 6  
63000 Cyberjaya, Malaysia

**Australia & New Zealand**

Ansell Limited  
Level 3, 678 Victoria Street  
Richmond, Vic, 3121  
Australia