

sempermed® 
derma plus

DERMA PLUS PRODUCED SAFETY

HIGH TEAR RESISTANCE – through double-immersion production process

PERFECT GRIP – through micro-rough surface

ANATOMICALLY PERFECT SHAPE – with broader back of the hand
and naturally curved finger



PRODUCED SAFETY

Manufactured while using the specially developed double-dipping process the Sempermed® Derma plus excels with good donning properties and an AQL of 0.65.

PRODUCT HIGHLIGHTS

- High tear resistance and reliable impermeability through double-immersion production process
- Anatomically perfect shape with broader and more comfortable back of the hand and naturally curved fingers
- Good grip through micro-rough surface texture, even in a moist environment



TECHNICAL DATA

sempermed® derma plus

TYPE	Sterile surgical glove for single use, lightly powdered			
MATERIAL	Natural rubber latex (NR)			
COLOUR	Natural white			
GLOVE SHAPE	Fully anatomical with rolled rim			
SIZE/OVERALL LENGTH as per EN 455-2	5½		min. 265 mm	
	6 and 6½		min. 270 mm	
	7, 7½ and 8		min. 280 mm	
	8½ and 9		min. 285 mm	
WALL THICKNESS measured in single layer in the palm area	0.19–0.24 mm			
BARRIER PERFORMANCE Tightness exceeding EN 455-1 (AQL 0.65 instead of AQL 1.5)	AQL 0.65			
FORCE AT BREAK exceeding EN 455-2 (≥ 12 N instead of ≥ 9 N)	≥ 12 N			
DURABILITY	5 years in original packaging when stored as per DIN 7716, ISO 2230			
STERILISATION	Radiation STERILE R			
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223			
PURPOSE as set out in MD Directive 93/42/EEC as set out in PPE Directive 89/686/EEC	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time) ¹⁾			
PACKAGING	Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical peel pack. in dispenser carton with sterilisation indicator: 50 pairs in transport carton with sterilisation indicator: 300 pairs			
ARTICLE NUMBERS	Size 5½	823152529	Size 7½	823152729
	Size 6	823152609	Size 8	823152809
	Size 6½	823152629	Size 8½	823152829
	Size 7	823152709	Size 9	823152909



CONTACT US!

Semperit Investments Asia Pte Ltd
8 Jurong Town Hall Road,
#29-03 to 06 JTC Summit, Singapore 609434
E-Mail: sempermed@semperitgroup.com
www.sempermed.com

Semperit Technische Produkte Gesellschaft m. b. H.
Division Sempermed
Austria · 1030 Vienna · Modecenterstrasse 22
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630
E-Mail: sempermed@semperitgroup.com · www.sempermed.com

¹⁾Please note the Sempermed Chemical Resistance list, available on www.sempermed.com/sempermed_informs/chemical_resistance

IMPORTANT NOTE: The latest product information is available under www.sempermed.com. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. Semperit does not accept any liability for incorrect use of the gloves. In case of doubt, obtain expert advice before use. The information and classification correspond to the latest status prior to printing. Technical details are average values from production and may vary in individual cases. Subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock.**